



Ric-Rac

CUSHION
by Martin Storey



YARN

Big Wool

5 x 100g balls

(photographed in White Hot 001)

NEEDLES

1 pair 10mm (no 000) (US 15) needles

Cable needle

EXTRAS

1 cushion pad approx. 57 cm (22½ in) square

TENSION

9 sts and 12.5 rows to 10 cm measured over st st, 13 sts and 12.5 rows to 10 cm measured over cable patt both using 10mm (US 15) needles.

SPECIAL ABBREVIATIONS

Cr4R = slip next st onto cable needle and leave at back of work, K3, then P1 from cable needle; **Cr4L** = slip next 3 sts onto cable needle and leave at front of work, P1, then K3 from cable needle; **Cr5R** = slip next 2 sts onto cable needle and leave at back of work, K3, then P2 from cable needle; **Cr5L** = slip next 3 sts onto cable needle and leave at front of work, P2, then K3 from cable needle; **C6B** = slip next 3 sts onto cable needle and leave at back of work, K3, then K3 from cable needle.

MB = make bobble as folls: (K1, P1, K1) all into next st, turn, P3, turn, K3, turn, P3, turn, sl 1, K2tog, pssso.

FINISHED SIZE

To fit a cushion pad approx. 57 cm (22½ in) square.